

AccuSway^{PLUS}

For Balance and Postural Sway Measurement

DESCRIPTION

AMTI's AccuSway^{PLUS} system is a complete solution for quantifying and evaluating human balance. AccuSway^{PLUS} was developed and built to be economical, portable, and easy to use.

The AccuSway^{PLUS}:

- provides a large 50 x 50 cm support surface
- can be used with a desktop or laptop computer
- plugs directly into the R232 port of the computer
- contains complete built-in electronics
- has an external trigger input that can be used to synchronize data
- weighs only 11.4 Kg
- is only 4.4 cm high

The AccuSway^{PLUS} measures the three forces (Fx, Fy, Fz) and the three moments (Mx, My, Mz)* involved in balance, providing outputs that allow easy computation of the center-of-pressure coordinates.

(*U.S. Patent #5,339,699)

A grid on the surface of the platform provides base-of-support coordinates which can be combined with the balance data, allowing the center-of-pressure (COP) to be plotted relative to the subject's position

SOFTWARE


AMTI's SWAYWIN95 postural sway measurement software acquires, analyzes, and plots data from the AccuSway^{PLUS} platform at a rate of 50 sets per second. SWAYWIN95 performs a comprehensive analysis of the sway data and presents many summarizing parameters. One particularly useful parameter, the area of a 95% confidence ellipse, may be used to track time related sway in balance due to the effects of drugs or therapy. Other results include the length of the path of the COP and the average velocity of the COP. Other parameters calculated by SWAYWIN95 are listed on the right.

SWAYWIN95 offers various plots of the sway data, including:

- color or black & white plots
- batch processing of plots
- the COP path on the force plate.
- actual boundaries of the 95% confidence ellipse.
- elliptical area vs. time plots during the test.
- Fast Fourier Transform (FFT) which shows a spectral analysis of the data.

SWAYWIN95 calculates statistical parameters using data from a single data file or from multiple sets of data files. SWAYWIN95's ability to calculate the minimum, maximum, average, and standard deviation across a set of data files allows the easy compilation of a baseline for any study.

SWAYWIN95 allows the user to export the COP and analysis data to comma delimited text files to be imported into spreadsheets like Excel, LOTUS, etc., for further analysis.



Parameters Calculated by SWAYWIN95

- dT Sampling interval
- N Number of data sets
- Xo, Yo Mean position of COP
- Xs, Ys Standard deviations of the coordinates of the COP
- Xm, Ym Mean deviations of the coordinates of the COP
- Rm Mean displacement of X, Y from Xo, Yo
- Rs Standard deviation of the displacement of X, Y from Xo, Yo
- Cc Correlation coefficient of X and Y
- Is Standard deviation of displacement along the axis of maximum standard deviation
- Iis Standard deviation of displacement along the axis of minimum standard deviation
- Theta Rotation required to bring the X-axis to the Is axis
- L Length of path of the COP
- VEL Average speed of COP along its path
- Ao Area included within the path of the COP
- A95 Area of the 95% confidence ellipse
- ABOS Area of the base-of-support

ISO 9001 CERTIFIED

AMTI

AGT Absolute Gauge Technologies

Presented by: Absolute Gauge Technologies

sales@absolutegauge.com; www.absolutegauge.com,



Toronto: 416 754 3168, Montreal: 514 695 5147, Toll Free: 1 888 754 7008

AccuSway^{PLUS}



For Balance and Postural Sway Measurement

AccuSway^{PLUS} Platform Specifications

Vertical Capacity:	1800 N centered on platform
Horizontal Capacity:	180 N centered
Mx, My Capacity:	140 Nm at corners
Height:	4.4 cm
Length x Width:	50 x 50 cm
Weight:	11.4 Kg
Data Rate:	50 Hz (each channel is read 50 times a second)
Interface:	RS232, 57.6K Baud, 8 bit, 1stop bit, no parity
Data Transmission:	Proprietary binary format
Power Supply:	12V DC, 110V or 220V input power supply included
Cable Length:	4.5 m
External Trigger:	Optoisolated signal input 0-0.9V (low) 2.5-15V (high), protected to ±18V
Interface Software:	AMTI's SWAYWIN95
Computer Requirements:	Desktop: Pentium 133 MHz processor, Windows95, 16 MB Ram Laptop: 233 MHz, 32 MB Ram



AccuSway^{PLUS} System Set-up



AMTI **ADVANCED MECHANICAL TECHNOLOGY, INC.**

Contents of this publication are subject to change without notice.