

NEW INTEGRATION! **MANUS™**

MR4 Now Supports MANUS Gloves for Enhanced Finger Kinematics

Seamlessly integrate anatomical data from the MANUS Gloves into your MR4 applications, enabling enhanced visualization and analysis of finger kinematics.

High-Precision Finger Tracking

- Record or livestream anatomical finger data in MR4 with unmatched accuracy.
- Capture 5 fingers, 15 joints, and 20 angles—per hand for in-depth motion analysis.

Modular & Scalable Design

- Top module designed for Ultium Motion sensors, ensuring seamless integration with Noraxon's ecosystem.

Effortless Usability

- Wear the gloves in seconds—no sensor placement hassle.
- Simple and efficient calibration procedure, allowing users to quickly set up and start capturing precise finger data.



Presented by: Absolute Gauge Technologies
sales@absolutegauge.com; www.absolutegauge.com,
Toronto: 416 754 3168, Montreal: 514 695 5147, Toll Free: 1 888 754 7008

NORAXON®

MANUS™

METAGLOVES PRO

Sensor sample rate	120Hz
Battery duration	Up to 3 hours (swappable)
Charging time	2 hours (USB-C)
Weight	154 grams
Wired communication	USB-C
Wireless communication	Bluetooth Low Energy 5
Wireless range	Up to 15 meters
Textile information	One size fits most (M/L)



 **AGT**

Presented by: Absolute Gauge Technologies
sales@absolutegauge.com; www.absolutegauge.com,
Toronto: 416 754 3168, Montreal: 514 695 5147, Toll Free: 1 888 754 7008

NORAXON®